

D.I WATER SYSTEM



FEATURES:

- Deionized water is superior to B.P. .specification..
 - No installation involved.
 - System operates on 240 Volts AC. (Optional-with pump)
 - Purified water immediately available, no storage required.
 - Convenient portable unit.
 - No regeneration required, simply replace resins refill.
- Optional: Pressure booster system

TECHNICAL SPECIFICATION:

Dimension : 660Wx720Hx610D (mm)
Weight: App. 30 kg
Water supply : Municipal water supply.
Flow rate : App. 120 litres per hour.

Model K8L Portable D.I system, built-in purity meter, pre-filter, post filter and stainless skid mount on caster wheels is design for used as a stand-alone system, suitable for up to two dialysis machines and for any other usage which requires 1 mega ohm quality de-ionized water